



**NEW!**

Validated under coordination  
of the Joint Research Centre



# See what you eat

DualChip® GMO Kit V2.0

**eppendorf**

# A new dimension in GMO screening

## DualChip® GMO Kit V2.0

The DualChip GMO Kit V2.0 is a multiplex screening tool for EU-authorized and non-authorized GMOs, successfully validated in collaboration with the Joint Research Centre (JRC) of the European Union\*.

The simultaneous detection of multiple GMO-specific elements, plant-species markers and control elements offers a time-sensitive and robust solution for GMO screening in food, feed and seed.

Optimized components assure best results and an incorporated software tool provides concise result reports with a few simple clicks.

### Screen and identify - the new content

Unique combination of plant-species, construct-specific and event-specific targets in one assay

Simultaneous detection of 30 transgenic elements for testing of GMOs in samples with unknown content

### Principle - as easy as this

- PCR-based amplification of target DNA with biotinylated primers
- Microarray-based detection of amplicons with specific DNA capture-probes
- Silver-based staining of hybridized targets using Eppendorf proprietary Silverquant® System
- Two independent tests per slide

### To get started - The "All-in-One" System


- Kit components include beside the microarrays all reagents necessary for amplification, labeling, hybridization including process controls
- Application-specific software tool is included, providing convenient analysis
- Comprehensive „All-in-One“ system for array processing, target detection and slide scanning.

**DNA extraction**  
from food or feed sample

**Amplification**  
via 3 different multiplex  
PCRs with biotinylated  
primers

**Hybridization**  
of the 3 PCR reactions  
on 1 DualChip GMO  
microarray

**Detection**  
with Silverquant®  
colorimetric detection



PDC 6	PHC	NHC	PDC 6	PDC 6
PDC 6	p35S	PDC 6	tNos	CaMV
PDC 13	NDC	pNos-nptII	PCRC S1	EPSPS-A
PDC 12	EPSPS-B	EPSPS-C	pat	NDC
PDC 11	PDC 6	cry1Ab-1	cry1Ab-2	cry1Ab-3
PDC 10	PCRC P1	bar	cry3Bb1	PDC 6
PDC 9	Rice	Sugar beet	Soybean	PCRC E1
PDC 8	Cotton	Maize	NDC	Potato
PDC 7	Brassicaceae sp.	PCRC S2	PDC 6	RRS
PDC 6	Bt11	PDC 6	Bt176	GA21
PDC 5	GT73	MON531	MON810	NDC
PDC 4	MON810	NDC	MON863	PCRC P2
PDC 3	MON1445	MON15985	NDC	C
PDC 2	PCRC E2	T45	PDC 6	NDC
PDC 1	PDC 6	NDC	NDC	NDC
PDC 6	PDC 6	PDC 6	PDC 6	PDC 6

**Figure 2:** Spotting pattern of probes and controls present on the DualChip GMO V2.0

\* <http://bgmo.jrc.ec.europa.eu/home/documents/report-JRC-EAT.pdf>




Figure 1: DualChip GMO Kit workflow

### Content

#### PCR S: Twelve screening elements:

p35S, tNos, pat, cry1Ab (3 gene varieties), cry3Bb1, EPSPS (3 gene varieties), bar, pNos-nptII

#### PCR P: Seven plant-specific markers for:

Maize, Soybean, Rapeseed, Cotton, Rice, Potato, Sugar beet

#### PCR E: Eleven event-specific markers for:

Bt11, Bt176, GA21, GT73, MON531, MON810, MON863, MON1445, MON15985, RRS, T45

#### Control elements:

CaMV virus contamination control  
 Positive PCR control (PCRC)  
 Positive and negative hybridization control (PHC, NHC)  
 Positive and negative detection control (PDC, NDC)

### Specifications

- Accuracy of 95% established in ring-trial study
- Proven sensitivity of 0.1% for all transgenic elements and 1% for all plant-specific markers
- All capture-probes (incl. controls) spotted in triplicate
- Tight monitoring of handling and results due to sophisticated control system for each step
- Standardized, sample-independent protocol
- High specificity from 2 independent sequence-specific steps
  - specific primers for the PCR reaction
  - specific capture probes for the hybridization
- Convenient software tool for concise report generation

DualChip <sup>®</sup> GMO Report			
Operator name: Silverquant detection Kit lot: DualChip-GMO Kit Box 1 lot: DualChip-GMO Kit Box 2 lot:	File name: X400027M X400119K X400228K	File name: 401817910008-1-Bt176.txt	Experiment date: 02.08.2008
Analysis date: 20.05.2008	Sample number: X400119K	Sample designation: 6.1% Bt176 CFM	Comments:
<b>GM event(s) present</b>			
<b>Genetic target elements</b>		<b>Specified GM events and hybrids based on event-specific target element</b>	
<b>Screening:</b>		<b>GM event</b>	<b>Event-specific</b>
pat	Detected	Bt176	Present
cry1Ab-1	Not detected	GA21	Absent
cry1Ab-2	Not detected	GT73	Absent
cry1Ab-3	Not detected	GT9 40-3-2	Absent
cry3Bb-1	Not detected	MON1445	Absent
EPSPS-A	Not detected	MON15985	Absent
EPSPS-B	Not detected	MON531	Absent
EPSPS-C	Not detected	MON810	Absent
RRS	Not detected	MON863	Absent
pat	Not detected	T45	Absent
Proo-nptII	Not detected	<b>GM hybrid</b>	<b>Event-specific</b>
Tnos	Not detected	CA21 x MON810	Absent
<b>Plant species:</b>		MON15985 x MON1445	Absent
Brassicaceae	Not detected	MON531 x MON1445	Absent
Cotton	Not detected	MON863 x MON810	Absent
Maize	Detected		
Potato	Not detected		
Rice	Not detected		
Soybean	Not detected		
Sugar beet	Not detected		
<b>Event-specific:</b>		<b>Specified GM events and hybrids based on pattern analysis</b>	
Bt11	Not detected	<b>GM event</b>	<b>Pattern analysis</b>
Bt176	Detected	CA51507	Absent
CA21	Not detected	DAS59122	Absent
GT73	Not detected	HT1	Absent
GT9 40-3-2	Not detected	NK903	Absent
MON1445	Not detected	T29	Absent
MON15985	Not detected	TP95S192	Absent
MON531	Not detected		
MON810 a	Not detected	<b>GM hybrid</b>	<b>Pattern analysis</b>
MON810 b	Not detected	CA51507 x NK903	Absent
MON863	Not detected	MON863 x NK903	Absent
T45	Not detected	MS1 x RF1	Absent
		MS1 x RF2	Absent
		MS6 x RF3	Absent
		NK903 x MON810	Absent
<b>Contamination control</b>		<b>Unspecified GM event</b> Detected screening target element(s)	
CaMV	Not detected		
<b>Process controls</b>			
Positive hybridization controls:	OK		
Negative hybridization controls:	OK		
Negative detection controls:	OK		
Positive detection controls:	OK		
PCR controls:	OK		
Array variability:	OK		

Figure 3: Example of DualChip GMO result table

# DualChip® GMO Kit V2.0

## Kit components

### Kit Box 1

- PCR control V2.0
- GMO primer mix S, P, E

### Kit Box 2

- 8 DualChip GMO screening slides V2.0 (16 individual tests)
- 16 Hybridization frames
- Aluminum sealing pads
- Genomic HybriBuffer
- NaOH Solution
- Hybridization control V2.0

### Kit Box 3

- Qiagen Multiplex PCR Master Mix for 60 reactions\*
- RNase-free H<sub>2</sub>O

DualChip GMO software tool for analysis of results and generation of ISO 17025-conform report accessible by download from: [www.eppendorf-biochip.com/downloads](http://www.eppendorf-biochip.com/downloads)

## Additional Hardware and Reagents

- Silverquant® Microarray Detection Kit for staining of 16 DualChip GMO microarrays
- Silverquant Microarray Scanning System for scanning and raw data analysis
- Thermomixer® comfort for simultaneous mixing and temperature control during hybridization
- Silverquant adapter for incubation of microarray slides during Silverquant staining

## Ordering information

Description	Order No. International	Order No. North America
DualChip® GMO Kit V2.0 (8 slides for 16 multiplex tests incl. reagents and software)	0038 004.032	956014200
Silverquant scanner system (software included)	6000 000.106	956030030
Silverquant detection kit	0038 005.000	956015073
Silverquant adapter	5368 000.037	956015090
Thermomixer comfort (without exchangeable thermoblock)	5355 000.011 (220V)	022670107 (120V) 022670158 (220V)
Exchangeable thermoblocks and accessories for 4 slides for hybridization experiments	5368 000.010	022670590

\* The QIAGEN Multiplex PCR Kit is intended for research use. No claim or representation is intended for its use to provide information for the diagnosis, prevention of treatment of a disease. This product is intended to be used for research purposes only.

The DualChip GMO Kit V2.0 is intended to be used for research purposes only.



Eppendorf Biochip Systems GmbH

22313 Hamburg · Germany · Tel: +49 40 538 01-501 · E-mail: [biochipsystems@eppendorf.de](mailto:biochipsystems@eppendorf.de) · [www.eppendorf-biochip.com](http://www.eppendorf-biochip.com)

**Application support** Europe, international: Tel: +49 1803 666 789 · E-mail: [biochipsupport@eppendorf.de](mailto:biochipsupport@eppendorf.de)

North America: Tel: +1 800 645 3050 ext. 2258 · E-mail: [support\\_na@eppendorf.com](mailto:support_na@eppendorf.com)

Asia Pacific: Tel: +60 3 8023 6869 · E-mail: [support\\_asiapacific@eppendorf.com](mailto:support_asiapacific@eppendorf.com)