pco.inspector high speed mobile digital high speed 10bit CMOS camera system with Tablet PC

- 636 fps at full resolution (1357 fps at VGA resolution)
- extremely mobile solution with Tablet PC LM1600
- high resolution (1280 x 1024 pixel)
- exposure time range 50 ns 5 s
- more than 2 h operation by battery
- image memory in camera (camRAM up to 4GB)
- easy to use tripod
- standard interfaces (IEEE1394, camera link)





pco.inspector high speed

This high speed 10bit CMOS camera inspection system comprises advanced CMOS and electronics technology. It is packed together with an excellent Tablet PC LM1600 into a case with batteries and tripod. Therefore, the inspection system can be used anywhere operated for more than 2 hours by its own batteries

With the new approach to integrate the image memory (camRAM) into the camera itself, it enables unmatched fast image recording with 1GB/s. The system features an excellent resolution (1280 x 1024pixel) and low noise. It consists of a compact camera with an external intelligent power supply. The image data are transferred via customer selectable standard data interfaces to a computer (IEEE 1394 ("firewire"), camera link). The available exposure times range from $1\mu s$ (50ns optional) to 5s. This digital CMOS camera system is perfectly suited for high speed inspection applications such as production line fault detection and motion analysis.

technical data: camera system

	unit	setpoint	pco.1200 hs
resolution (hor x ver) ¹	pixel		1280x1024
pixel size (hor x ver)	μm ²		12.0 x 12.0
sensor format/	mm²/		15.36x12.29/
diagonal	mm		19.67
peak quantum	%	@ 520nm	27
efficiency		typical	
full well capacity	e-		63 000
image sensor			MT9M413
dynamic range	dB	@ CMOS	59.6
		camera	
dynamic range A/D ²	bit		10
readout noise	e ⁻ rms	@ 66MHz	85
imaging frequency,	fps	@full frame/	636/
frame rate		@ROI VGA	1357
pixel scan rate	MHz	dual speed	66 / 86
A/D conversion factor	e ⁻ /count		30
spectral range	nm		2901100
exposure time	S		1µs5s
			(50ns5s opt.)
anti-blooming factor		typical	no blooming
smear	%		no smear
binning horizontal	pixel		1
binning vertical	pixel		1
dark current	e ⁻ /pixel·s	@25 °C typical	5900
region of interest	pixel	horizontal	steps of 10
		vertical	steps of 1
interframing time	ns	@FWHM ³ and	70
(PIV mode)		100% fullwell	
		signal	



technical data: camera system

non linearity	%	full temperature range	<2
uniformity darkness DNSU ⁴	e ⁻ rms	@ 90% center zone	<700
uniformity brightness PRNU ⁵	% rms	typical	0.6
trigger, auxiliary signals		internal/ external	software / TTL level
power consumption	W	typical/ maximum	25 / 40
power supply	VAC		90260
mechanical dimensions camera (w x h x l)	mm ³		84 x 66 x 175
mechanical dimensions case (w x h x l)	mm ³		178 x 425 x 500
weight	kg	complete system	12.9
operating temperature range	°C		+5+40
operating humidity range	%		1090
storage temperature range	°C		-20+70
optical input			Nikon f-mount, c-mount
data interface			IEEE1394, camera link, ethernet
CE certified			yes
			link, ethernet

software:

Camware software for camera control, image acquisition and archiving of images in various file formats, WindowsXP and later, 32bit-dynamic link library (DLL) is available for user customisation and integration on PC platforms (software development kit - SDK), software is operational in either single mode or with builtin recorder functions, drivers for popular third party software packages are available (see website)

options: CMOS image sensor in color version

custom-made versions

camRAM available in: 1GB, 2GB and 4GB

trolley for simple transport and as solid base for operation

- horizontal versus vertical Analog-to-Digital-converter full width half maximum
- dark signal non-uniformity
- photo reponse non-uniformity



frame rate table [frames per second]

pixelclock	66MHz	86MHz
exposure time	1/fps / <1/fps	1/fps / <1/fps
1280x1024 pixel (full frame)	488 / 486	636 / 634
1280x512 pixel	977 / 969	1272 / 1263
1280x256 pixel	1953 / 1923	2545 / 2506
1280x128 pixel	3906 / 3788	5090 / 4936
1280x64 pixel	7813 / 7353	10180 / 9581
1280x32 pixel	15625 / 13889	20360 / 18098
1280x16 pixel	31250 / 25000	40720 / 32576

(Regions of Interest are free selectable, this table just shows the relation between vertical resolution and image rate)

Tablet PC - Motion Computing LM1600



technical data: Tablet PC LM1600

processor		@ 1.5 GHz	Intel [®] Pentium [®] M processor LV758
operating system			MS [®] Windows [®] XP Tablet PC Edition
RAM	GB		1
harddisk			60
display		XGA TFT LCD View Anywhere®	12"
weight	kg		1.42
mechanical dimensions Tablet PC (w x h x l)	mm ³		296 x 240 x 18.7-22
operating temperature range	°C		+5+35
operating humidity range	%		2080
storage temperature range	°C		-20+65
CE certified			yes

Further technical information can be found on: www.motioncomputing.com or www.thetabletstore.de



View into a bottling machine, which fills beer bottles. The uptake of the empty bottles can be seen in the background (recorded at 636fps).

...with friendly permission of G. Schneider & Sohn GmbH, Kelheim, Germany



areas of application

■ high speed inspection ■ motion analysis in sports ■ fault detection in high throughput production lines ■ motion analyis of animals ■ material testing ■ tensile tests ■ airbag inflation ■ fast flow visualisation ■ spray analysis ■ hydrodynamics ■ fuel injection ■ fast events in nature and medicine

contact

PCO AG Donaupark 11 93309 Kelheim, Germany

fon +49 (0)9441 2005 50 fax +49 (0)9441 2005 20 info@pco.de www.pco.de

pco.inspector highspeed 12/2005 - preliminary subject to changes without prior notice © PCO AG, Kelheim

The Cooke Corporation 6930 Metrolex Drive Romulus, MI 48174 USA

fon +1 248 276 8820 fax +1 248 276 8825 info@cookecorp.com www.cookecorp.com

