

ExpellPlus™

Superior performance tips

Due to special polymers and an extremely smooth inside surface, the ExpellPlus™ tips from Capp® virtually eliminate sample retention and reduce the binding of DNA, RNA and Proteins to a very minimum. Perfectly matched with the CappAero™ and CappTronic™ pipettes, they provide a superior pipetting performance, which means savings of valuable samples and improved pipetting accuracy and precision.

Plus™ low retention properties

Tests show that the ExpellPlus™ tips deliver volumes within 0.1 % of the set volume, versus 0.5-0.7 % for standard tips. This significant difference in sample retention also is clearly visible when pipetting a dyed liquid (see below).



Standard Tip



ExpellPlus™ low retention

The Expell™ tips without the Plus™ low retention properties are of the same high manufacturing standards, maintaining strict tolerances with ideal wetting properties and high transparency and purity levels.

Cost saving ReLoad system

A unique ReLoad system makes re-filling of racks easier than ever, and is designed to save cost, and waste as well.



BluCapp™ 15 & 50 ml Centrifuge Tubes

Designed to cover the everyday need for sample storage and centrifugation in the lab, the BluCapp™ feature easy to read graduations as well as a large writing area (20x50mm). The leak proof caps have large ridges to provide the optimal grip when wearing gloves.

Combined with the ¾ turn closure, it significantly increases speed and allows one-handed operation. Made of polypropylene, the BluCapp™ centrifuge tubes are autoclavable, and are tested to endure up to 6000 RCF (15ml) and 9400 RCF (50ml).



ExpellPlus™ Filter tips

The ExpellPlus™ tips also comprise of a very affordable and high quality range of filter tips, the use of which eliminates any risk of cross contamination. The ExpellPlus™ filter tips are tested and certified to conform to strict bio-load specifications:

- No detectable nucleic acid contamination
- No detectable PCR inhibitors
- No detectable endonucleases
- No detectable endotoxin
- Sterile




Easily accessible, space and waste saving

The racked Expell™ tips are supplied in handy dispenser boxes occupying minimal space in the bench, with easy access to the next rack to be used. Due to innovative packaging solutions, the racked tips take up less storage space than conventional tips and are waste saving at the same time.




We take pipetting to perfection.

Expell™ and ExpellPlus™ Filter Tips



Cat. No.	Description	Qty	Packaging
5030030	ExpellPlus™ 10µl FT	10 racks of 96	Racked/Sterile
5030060	ExpellPlus™ 10µl FT, long	10 racks of 96	Racked/Sterile
5030062	ExpellPlus™ 20µl FT	10 racks of 96	Racked/Sterile
5030066	ExpellPlus™ 100µl FT	10 racks of 96	Racked/Sterile
5030090	ExpellPlus™ 200µl FT	10 racks of 96	Racked/Sterile
5130150	Expell™ 1000µl FT	8 racks of 60	Racked/Sterile


Expell™ and ExpellPlus™ Pipette Tips



Cat. No.	Description	Qty	Packaging
4030005	ExpellPlus™ 0.2-60µl*	10 racks of 384	Racked/Sterile
4030006	ExpellPlus™ 0.2-30µl*	10 racks of 384	Racked/Sterile
<i>* special for use with 16-channel pipettes from Capp.</i>			
5130010	Expell™ 10µl	1000	Bulk
5030010	ExpellPlus™ 10µl	1000	Bulk
5030015	ExpellPlus™ 10µl	10 racks of 96	Racked
5030020	ExpellPlus™ 10µl	10 racks of 96	Racked/Sterile
5030035	ExpellPlus™ 10µl Reload	20 cards of 96	Reload/Sterile
5130040	Expell™ 10µl, long	1000	Bulk
5130045	Expell™ 10µl, long	10 racks of 96	Racked
5130050	Expell™ 10µl, long	10 racks of 96	Racked/Sterile
4130075	Expell™ 200µl, yellow	1000	Bulk
5130070	Expell™ 200µl	1000	Bulk
5030070	ExpellPlus™ 200µl	1000	Bulk
5030075	ExpellPlus™ 200µl	10 racks of 96	Racked
5030080	ExpellPlus™ 200µl	10 racks of 96	Racked/Sterile
5030095	ExpellPlus™ 200µl	10 cards of 96	Reload/Sterile
5130100	Expell™ 300µl	1000	Bulk
5130105	Expell™ 300µl	10 racks of 96	Racked
5130110	Expell™ 300µl	10 racks of 96	Racked/Sterile
4130135	Expell™ 1000µl, blue	1000	Bulk
5130130	Expell™ 1000µl	1000	Bulk
5130135	Expell™ 1000µl	16 racks of 60	Racked
5130140	Expell™ 1000µl	16 racks of 60	Racked/Sterile
5130160	Expell™ 5000µl	250	Bulk

BluCapp™ Centrifuge Tubes

Cat. No.	Description	Qty	Packaging
5100050	BluCapp™ 50ml	50	Bulk/Sterile
5100015	BluCapp™ 15ml	50	Bulk/Sterile



Capp A/S
 Landbrugsvej 10 • DK 5260 Odense S
 Denmark
 tel. +45 6613 6140 • fax. +45 6613 2770
www.cappbrand.com • info@capp.dk

(distributor)